Measurement of β -ray spectra

Ana C. Fabela Hinojosa

School of Physics and Astronomy, Monash University

Experiment performed: Tuesday 18th August, 2020

Abstract

Using a thin lens magnetic spectrometer, we measure the momentum spectrum of electrons emitted as β^- rays from a radioactive source of $^{137}\mathrm{Cs}$. The detected momentum of the radiated electrons is defined by the spectrometer's adjustable magnetic lens current and k a proportionality constant dependent on the geometry of the apparatus. The magnetic field of the lens is varied by changing the current passing though the lens coil which has the effect of modifying the trajectories of the electrons, focusing electrons with specific momenta onto the detector allowing us to measure their intensity. By converting the measured momentum to energy we are able to fit our data to a linear model based on the Fermi-Kurie plot. We find that the value of the kinetic energy of the nuclear transition is $T=0.520\pm0.044$ MeV which is in agreement with the accepted value of T=0.512 MeV[1].

1 Introduction

When Henri Becquerel first observed β —radiation, he determined that the observed radiated particle satisfied the same mass-to-charge ratio as the electron, the particle that J.J Thompson had discovered in 1897.

Later experimental results showed that electrons observed as β —rays are emitted with a continuous range of kinetic energies up to a maximum value[2]. The discovery of a continuous distribution rather than a discrete predictable value for the electron kinetic energies took physicists by surprise since it violates energy, momentum and angular momentum conservation laws.

It was in fact this counterintuitive result which led Wolfgang Pauli to propose in 1930 that the observed anomalies must be due emission of a yet unknown particle. In 1934 Enrico Fermi called this massless and undetectable particle a "neutrino" and developed an advanced theory of beta decay based on it, but the neutrino was not experimentally observed until much later in 1956.[3] Due to symmetry considerations, the lepton emmited together with the electron during β^- decay is called an antineutrino.

The process we currently know as β^- decay describes a neutron in a parent nucleus desintegrat-

ing into a proton in a daughter nucleus, an electron and an antineutrino.

In a β -event, both nuclides (nuclear species) have the same number of nucleons. This means that the daughter nucleus will not experience a substantial change in kinetic energy (recoil) due to the decay event. Leaving most of the desintegration energy available to be carried-off by the leptons as kinetic energy.

A parent nucleus has a given initial energy w. The avalilable kinetic energy of the system is equal to the decrease in mass energy due to the creation of the radiated leptons:

$$T = w - mc^2, (1)$$

where m is the difference in mass between the daughter and parent nuclides.

The observable count of β^- electrons as a function of energy is now described by the Kurie–Fermi Theory of β^- decay. The theory defines the distribution of the available nuclear desintegration energy to be carried-off as kinetic energy by both leptons. The observed maximum value in the Kurie–Fermi plot represents the count of electrons that take carry-off the maximum possible kinetic energy whilst the antineutrinos carry close to zero kinetic energy from the transition.

2 β decay theory

In this experiment we measure the momentum spectrum of emitted β -rays from a radioactive source of ¹³⁷Cs. ¹³⁷Cs decays to an excited state of ¹³⁷Ba. This transition occurs with a probability of 94.6%[1]. The maximum energy value of this decay event is known to be T=0.512 MeV.

Our experimental apparatus is a thin magnetic lens spectrometer. The operation of β spectrometers depends on the behaviour of electrons sub-

ject to magnetic fields. The magnetic field of the spectrometer lens is varied by changing the current passing though the lens coils which has the effect of modifying the trajectories of the electrons by the Lorentz force,

$$\vec{F} = e^- \vec{v} \times \vec{B} \tag{2}$$

The result is a set of electrons focused onto the spectrometer detector all with a specific range momenta (other electrons in the momentum spectrum undergo chromatic aberration)

References

- [1] Monash SPA. 4.4 measurement of β -ray spectra, 2020.
- [2] C A Moyer R A Serway, C J Moses. Modern Physics, third edition. 2005.
- [3] Carl R Nave. Beta radioactivity, 2001.